
## T327
## COMBI-BOX™

Made of pvc

The Simport® Combi-Box™ can be used not only as a storage rack but also as a workstation. The white base will accept all 96- and 384-well plates and an easy to remove transparent cover allows easy viewing of contents. Being only 40 mm high, the stackable Combi-Box™ saves space on the lab bench and on refrigerator or freezer shelves. For single tubes as well as strips, use the Combi-Rack™ (T327-1) which can hold up to 96 tubes or 12 strips of 8.

| Cat. # | Color | Qty/Cs |
|---|---|---|
| T327 | White base | 5 |



- Plate can easily be inserted and removed
- Alphanumeric identification of each position
- Can accept single tubes, strips and plates up to 96 wells
- Stackable

## T327-1
## COMBI-RACK™

Made of polypropylene

The Simport® Combi-Rack™ is an innovative support that can hold up to 96 PCR tubes or 12 strips of 8 tubes with caps. Each hole is identified with an alphanumeric numbering system for identifying tubes. The grid stands on 4 legs and can be placed on a lab counter or in a refrigerator or freezer shelf. Made of polypropylene, it can easily withstand temperatures from -80 °C to +121 °C. It is also ideal for carrying and storing, freezing and transporting reagents and specimens.

For storage, simply place the Combi-Rack™ in the T327 Combi-Box™ and place cover.

| Cat. # | Color | Qty/Cs |
|---|---|---|
| T327-1 | Blue | 5 |

## T328-96
## PCRack™

Made of polypropylene

This convenient space saving rack was designed especially for storing and working with PCR samples. The PCRack™ will accept all models of 0.2 ml tubes, along with strips of 8 or 12 tubes. 96-well PCR plates can also be accommodated.

The PCRack™ can be horizontally attached to each other in order to build-up any configuration you desire. With the cover on, they are easily stackable one on top of another. Thanks to these special features, efficiency is highly improved allowing you to carry a multitude of tubes and/or strips at the same time.

| Cat. # | Color | Qty/Cs | Cat. # | Color | Qty/Cs |
|---|---|---|---|---|---|
| T328-96B | Blue | 20 | T328-96R | Red | 20 |
| T328-96G | Green | 20 | T328-96Y | Yellow | 20 |
| T328-96O | Orange | 20 | T328-96AS | Assorted* | 20 |
| T328-96P | PinK | 20 | | | |

*Blue, green, orange, pink, yellow.



## S500-80
## UniRack™

Made of polypropylene

The UniRack™ offers the laboratory a support far more versatile and easy to use than any other rack available today. It is designed to use minimum counter space while offering maximum flexibility. Made of polypropylene, it allows great resistance to various chemicals used in laboratories. For further information, see page 142.

| Cat. # | Color | Qty/Cs | Cat. # | Color | Qty/Cs |
|---|---|---|---|---|---|
| S500-80B | Blue | 10 | S500-80R | Red | 10 |
| S500-80G | Green | 10 | S500-80Y | Yellow | 10 |
| S500-80O | Orange | 10 | S500-80AS | Assorted* | 10 |
| S500-80P | PinK | 10 | | | |

*Blue, green, orange, pink, yellow.

| Lid Cat. # | Color | Qty/Cs |
|---|---|---|
| S501-80 | Transparent | 10 |

